A decorative graphic consisting of a trail of small, light blue particles that curves across the top half of the slide, resembling a comet tail or a molecular path.

Enzymatic Oligonucleotide Synthesis Process Flow and Substance Impurity Profile

May 14th 2024

Catalyzing Innovation Through Engineered Enzymes

Powerful Underlying Science

CodeEvolver® directed evolution platform underpins exquisite enzyme engineering capabilities



Enzymes from Nature



Deep Technical Expertise

20+ years evolving enzymes supporting **small molecule API production and life science applications**



Highly differentiated products



Custom Solutions for Biocatalysis

Over 50 commercialized enzymes supporting small molecule API production



Pharmaceutical Manufacturing




RNA Manufacturing

Potential Game-Changing Technology

ECO Synthesis™ manufacturing platform poised to disrupt the RNAi synthesis market


Enzymatic Solutions for siRNA Synthesis are Increasing

Chemical




Traditional siRNA Synthesis using Phosphoramidite Chemistry (PAC)

Hybrid



PAC fragments + Ligation

Enzymatic



ECO Synthesis™ (Enzyme-Catalyzed Oligonucleotide Synthesis) Manufacturing Platform

ECO Synthesis™ Manufacturing Platform: Process Flow and Impurity Profile

Spotlight Presentation Luncheon 1

Derek Gauntlett, MBA, Director, ECO Process Development

Two Enzymatic Approaches for Large-scale siRNA Synthesis

TIDES Talks in the Exhibit Hall

Thursday, 16 May 2024 10:10 - 10:20 AM

Mathew Miller, PhD, Director, Life Sciences & RNA Technology

A decorative graphic consisting of a trail of small, light blue particles that curves across the top half of the slide, resembling a comet tail or a molecular path.

ECO Synthesis™ Manufacturing Platform

Derek Gauntlett

CODEXIS®

ECO Synthesis™
Manufacturing
Platform:
Modularity for
Maximum
Versatility

Core Technology:
Sequential
Incorporation
of Modified
Building Blocks

Reagents:
NQPs & Starter
Oligo

Conjugation:
Incorporation of
Targeting Moieties

Ligation:
(dsRNA Ligase
Program)

ECO Synthesis™ Manufacturing Platform – Key Elements




Starter Oligo

- Substrate from which synthesis initiates
- Synthesized via phosphoramidite chemistry (PAC) or enzymatically




Nucleotide Quad Phosphate (NQP)

- Modified nucleotide monomers with blocking group
 - Phosphoramidite analog
- Synthesized chemically or enzymatically



Extension Enzyme


- Engineered to incorporate 2'-modified NQPs



Deblocking Enzyme

- Engineered to remove blocking group

ECO Synthesis™ Manufacturing Platform - Overview




Key Performance Indicators for ECO enzymes:

- Non-native Substrates
 - Engineered for modified nucleotides
- Promiscuity
 - Minimal bias on 300k+ possible substrates
- Productivity
 - Coupling efficiency ~99% in hours
 - Minimal substrate & product inhibition
- Robustness
 - Stability
 - Manufacturability


Targeting Moieties for Therapeutic Relevance

ECO Synthesis™ Manufacturing Platform Core Cycle




Enzymatic conjugation of targeting moieties

4 Conjugate




Targeting Moieties for Therapeutic Relevance



Enzymatic conjugation of targeting moieties

4 Conjugate



ECO Synthesis™ manufacturing platform enables quantitative incorporation of modified nucleotides with reactive and sterically demanding bulky side chains at 5', 3' and intra oligonucleotide.

Chemical Synthesis vs. ECO Synthesis™ Manufacturing Platform

Phosphoramidite Chemistry

Limited Scalability

- Best suited for bench top scale due to limited, single-digit kilogram batch sizes
- Capacity will be challenged to support future demand

Toxic Solvent Use

- Requires large volumes of inorganic solvent (acetonitrile) with high disposal costs
- Likely future supply chain limitations & price volatility

Low Purity

- Inefficient for longer RNAs
- Significant impurities from complex protection / deprotections

High Cost

- High-cost infrastructure investment & raw materials
- High purification costs
- Expensive waste disposal



ECO Synthesis™ Manufacturing Platform

Scalable

- ✓ Enzymes and flow-based process
- ✓ Dual reactor setup enables larger batch sizes compared to PAC

Reduced Waste

- ✓ Aqueous reactions
- ✓ Significantly decreases chemical waste streams
- ✓ Path to enzymatically created monomers

High Purity

- ✓ Higher purity driven by milder conditions during synthesis

Valuable Economics

- ✓ De-bottlenecks current supply constraints with increased scale & efficiency
- ✓ Integrates with existing manufacturing facilities that have purification capabilities

Progress Achieved on the ECO Synthesis™ Manufacturing Platform

Demonstrating the Power of Combing Enzyme Engineering & Process Development in One Year

Incorporation efficiency: **~88%**

FLP purity at N+4 (max. 4) = **~60%**

Accepting mod. nucleotides:
2'-OMe, 2'-F

Incorporation efficiency: **~92%**

FLP purity at N+4 (max. 6) = **72%**

Accepting mod. nucleotides:
2'-OMe, 2'-F, PS

Incorporation efficiency: **>98%**

FLP purity at N+4 (max. 14) = **92.7%**

Accepting mod. nucleotides: **2'-OMe, 2'-F, PS, 2'-H, "others"**

Conjugation chemistry

May 2023
TIDES USA

Nov 2023
TIDES EU


Dec 2023

May 2024
TIDES USA

Gram-scale synthesis
(N+6 w/ fully modified
nucleotides)




ECO Synthesis™ Current Process Overview: 4-Step Unit Operation



- Immobilized polymerase to elongate the substrate in a controlled manner.
 - Fixed packed bed reactor with co-immobilized enzymes on resin
 - Recirculation flow model
 - One reservoir system
- Immobilized phosphatase to remove blocking group from substrate and NQP.
 - Fixed packed bed reactor with immobilized enzyme on resin.
 - Linear flow model
 - Two reservoir system
- Ultra filtration diafiltration (UFDF) to increase concentration and desalt
 - Desalt only performed as needed - compound dependent.
 - TMP Flux model
 - Cassette style system
- Enzyme to add conjugation moieties in batch
 - Batch reactor system

ECO Synthesis™ Liquid Flow Cycle – Current Research Scale



Key Process Conditions

1. Scale: 2mM (Starter) / 1.5 e.q. NQP
2. Volume Start: 9 mLs

3. Extension Residence Time: 0.10 hours
4. Extension Column Size: 3g resin / column

5. Deblock Column Size: 3g resin / column
6. Deblock Residence Time: 0.10 hours


Note: UFDF performed using 1kDa MWCO cassette at 5th and 13th extensions

Process conditions are not finalized and are under optimization


ECO Synthesis™ - Process by the Numbers

- Two fully 2'-modified RNA oligonucleotides were synthesized using ECO Synthesis™ manufacturing platform demonstrating the technology's capability to synthesize different sequences.

Compound	Number of Extensions	RNA Starter Length	Added RNA Base Types	Avg. Coupling Efficiency (%)	Avg. Coupling Time (H)	Avg. Deblock Time (H)
Compound 1 w/ 3' Alkyne	14	7	2'F / 2'OMe	98	11.1	2.5
Compound 2	11	7	2'F / 2'OMe	98	6.0	2.5



ECO Synthesis™ Impurity Profile MS - Compound 1 w/ 3' Alkyne



Key Impurity Takeaways

-Notable PAC impurities not present by MS

- Abasic, acetyl, IBU, cNet, and methylation

-Oxidation of PS->PO is the predominate impurity generated using ECO Synthesis™ today.



-Overextension products are second largest impurity present.

-Enzyme / Process fixes are currently underway for both sets of impurities.


Figure. TIC of Compound 1 w/ 3'Alkyne

ECO Synthesis™ Process Design – Commercial Scale


Extension / Conjugation



Deblock



Desalt



Envisioned commercial scale utilizes automation and throughput for improved volumetric productivity.

To Purification

Forward Looking for ECO Synthesis™ Manufacturing Platform


Codexis is Positioned for Significant ECO Synthesis™ Improvements in 2024

- ECO Synthesis™ manufacturing platform has made significant progress since TIDES USA in May 2023
- Continue process development
 - Holistic process design space created using a QbD approach to push >99% conversion efficiencies
 - Phosphorothioate oxidation control under sequential ECO Synthesis™ platform
 - Decreased coupling times for challenging sequences for volumetric productivity improvements
 - Pushing research scale towards fully automated commercial scale processes
 - Evaluate therapeutically relevant conjugate moieties
- Continue enzyme engineering

Codexis Center of Excellence for Enzymatic RNA Synthesis


Your Pathway to Lower Cost and More Scalable RNA Manufacturing

Hybrid



RNA Ligase Services

Enzymatic



ECO Synthesis™
Manufacturing Platform



Booth #628