

# KNOW LABS


## TRANSFORMING NON-INVASIVE MEDICAL DIAGNOSTICS

**ATTD Tech Fair**

March 7, 2024

# Know Labs

We are building a  
**non-invasive**  
**continuous glucose**  
monitoring device.





- Dimensions (w x d x h) 47 x 76 x 20 mm
- Weight. 100 g

RF Dielectric Spectroscopy

No Needles

100% Non-invasive

Continuous Monitoring

Real-time Data

No-Consumables

Rechargeable

Wearable


Adhesive or Strap

AI/ML-Powered Algorithms

100+ Potential Applications

# RF Dielectric Sensor

IP-PROTECTED | INCLUDED IN THE KNOWU



Energy Field Dimensions at 2,500 MHz

- w x d x h: 2.54 x 2.00 x 1.27 cm

- 5.4 cm<sup>3</sup> of volumetric data

**Antenna Array** that emits and captures radio wave signals in the microwave spectrum and generates an “**Energy Field**”, collecting “volumetric data”


**RF Generator** enables frequency sweeps from 300 to 4,400 MHz, at various intervals, **1.5M data points collected per hour = >400 per second**

**6 Key Parameters, customizable with each sweep:** power, frequency range, frequency step, dwell time, and antenna permutations = >30,000 combinations



# Overcoming Limitations

Our RF Dielectric Sensor sweeps the spectrum to collect **volumetric high-resolution data at high speed**, overcoming a limitation faced by fixed wavelength optical sensors.

## Optical Sensors (400nm – 1000nm)






-- All Health Wearable Devices --



# High-Resolution Data for Multiple Analytes...


*In Vitro* Experiment

Know Labs Sensor Proprietary Spectral Signal  
PBS Solutions with 200 mg/dL analyte concentration




# .. and for *In Vitro* and *In Vivo* Applications

## *In Vitro* Glucose Solutions Readings



**IN VITRO:** ADC Voltage (y-axis) measuring voltage variance based on glucose concentration and frequency sweeps

## *In Vivo* Glucose Readings Over 3 Hour Test



**IN VIVO:** ADC Voltage (y-axis) measuring voltage variance based on dielectric permittivities of blood glucose and frequency sweeps

# Validation: Stability, Repeatability and Accuracy

	2021		2022		2023		...
Manuscript	Proof of Principle with Mayo Clinic	Exploratory Clinical Study	Proof of Concept Clinical Study	Technical Feasibility Study	New Algorithm Refinement Study	Data Preprocessing Techniques Study	
Description	Demonstrated the accuracy of Bio-RFID sensor in quantifying different analytes <i>in vitro</i> (liquid solution).	First indication that Bio-RFID could be an accurate alternative to FDA-cleared glucose devices.	Proof of concept ability to quantify blood glucose non-invasively using RF.	Demonstrates Bio-RFID can deliver stable, repeatable results in measuring blood glucose levels.	Algorithm refinement in the non-invasive detection of blood glucose using Bio-RFID technology.	Improvement in machine learning model accuracy on an expanded mixed cohort dataset.	
Accuracy	100% <i>in vitro</i> accuracy	MARD <u>5.3%-6.7%</u>	MARD 19.3%	MARD 20.6%	MARD 12.9%	MARD <u>11.3%</u>	
# Participants	na	2	1	5	5	13	
# Datasets	na	3	22	106	106	366	
# Bio-RFID datapoints	na	<u>1.5M</u>	~183M	~430M	~430M	<u>~1.7B</u>	
# Reference Observations	na	75	~383	~1,555	~1,555	~3,311	




# ATTD 2024 Poster: ~11% MARD in Normal & Hyper

## Non-Invasive Blood Glucose Monitoring in People with Diabetes Using an RF Sensor and Venous Blood Comparator

D. Klyve, J. Anderson, K. Currie, C. Ward, K. Pandya, V. Somers





- 30 participants with [prediabetes and Type 2 diabetes](#)
- [Venous blood](#) as a comparative reference
- 3-hour [Glucose Tolerance Test](#) (GTT)

Glucose Range (mg/dL)	n	MARD (%)	±15%	±20%
Hypoglycemic (<70)	4	9.5 ± 8.3	75.0 ± 4.2	100.0 ± 0.0
Normoglycemic (70-180)	99	11.0 ± 2.7	75.8 ± 0.8	83.8 ± 0.7
Hyperglycemic (>180)	27	11.5 ± 3.1	66.7 ± 1.8	85.2 ± 1.3
<b>Total</b>	<b>130</b>	<b>11.1 ± 2.1</b>	<b>73.8 ± 0.8</b>	<b>84.6 ± 0.6</b>



100% of estimations in Risk Grades  
A and B (82.3% in A, 17.7% in B)

# KnowU Development Roadmap

Proof-of-Concept	Generation 0	Generation 1	KnowU
			
<ul style="list-style-type: none"><li>• <b>Exploratory</b> design</li><li>• Multiple components wired to each other</li><li>• Signal testing purpose</li><li>• 2 x 3 ft board</li></ul>	<ul style="list-style-type: none"><li>• <b>Miniaturized</b> format</li><li>• Wired connection to power source and data capture</li><li>• Restricted to <b>laboratory controlled environment</b></li></ul>	<ul style="list-style-type: none"><li>• <b>On-the-go</b> form factor</li><li>• Place your palm or arm for an on-demand, non-invasive blood glucose level</li><li>• Computer mouse size</li></ul>	<ul style="list-style-type: none"><li>• <b>Wearable</b> form factor</li><li>• Continuous monitoring</li><li>• 86% smaller and 68% lighter than previous Generation</li></ul>

# Generation 1 Details: Research Form Factor



## Gen 1

- Computer mouse size
- Research Lab in your pocket to accelerate data collection



Using Gen 1 Prototype on Forearm











Using Gen 1 Prototype on Hand

# KnowU Details: Wearable Form Factor



## KnowU Product Specifications

-  **Wearable**  
Adhesive, or strap for wrist or arm. Clip system allows for easy sensor attachment.
-  **Battery**  
24-hour continuous use. Rechargeable.
-  **Algorithm**  
Multiple algorithms in development. See Know Labs Scientific Validation mobile app page.
-  **Mobile app**  
Designed to companion the KnowU CGM device. Enables advanced research.
-  **Data Storage**  
256MB flash storage, 40 hours of offline research data.
-  **Connectivity**  
BLE, WIFI, USB-C.
-  **Processor**  
Ultra-low power consumption, onboard computer. Wi-Fi & BLE capabilities, built-in ML.
-  **Firmware**  
Custom firmware supports the entire tech stack.



KnowU vs. Gen 1

86% smaller

68% lighter

# KnowU Real-life Application



Worn with Strap



Worn with Adhesive

Mobile App



# 2024 & Beyond

## Clinical Data Collection

More participants with diverse backgrounds, including people with T1D, T2D, and pre-diabetes.

## Algorithm Refinement

More data equals better accuracy – potential to accelerate time to market with calibration.

## Hardware Development

Device adjustments and/or further miniaturization

## Scientific Validation

External research institutes to further validate our technology and support FDA application.

## Intellectual Property


#1 worldwide in non-invasive blood glucose monitoring (>300 patents filled and pending)

## Increase the generalizability of the Bio-RFID


Must be safe and deliver required level of accuracy under any condition and regardless of the patient, as determined by its intended use

# Visit us at Booth #30


Try the KnowU  
and Scan your  
RF Signature



Air  
Seattle-WA



Air  
Florence-ITA



Body  
Leo Trautwein



# THANK YOU

Leo Trautwein  
Chief Commercial Officer  
[leo@knowlabs.co](mailto:leo@knowlabs.co)

